



Prodent Ltd

Quality Dental & Surgical Instruments, Imported from Pakistan, Trusted in Mauritius

Dental Instrument — Dental mirror

Category product list

Dental mirror



Anti fog mirror with handle

Anti fog mirror Front surface mirror Stainless steel Rust free Removeable handler



Dental mirror 12 piece

Professional dental mirror designed to provide a clear view of the oral cavity during dental examinations and procedures. The mirror offers clear reflection and comfortable handling for everyday clinical use. Key Features:- Stainless steel. Rust-free. Magnet-free. Clear reflection. Suitable for professional dental use. Easy to clean and sterilize.



Micro Mirror (4 Sizes)

The Dental Micro Mirror is a precision diagnostic and endodontic instrument designed to provide exceptional visibility in hard-to-reach areas of the oral cavity. Featuring a high-reflective mirror surface and a lightweight ergonomic handle, it delivers clear indirect vision for detailed examination and precision treatment. Available in 3 mm, 5 mm, 2×7 mm, and 3×9 mm mirror sizes, this instrument is ideal for endodontic procedures, restorative dentistry, and microscopic clinical applications where enhanced visibility is essential. Available Sizes : 3 mm – Round Micro Mirror 5 mm – Round Micro Mirror 2×7 mm – Rectangular Micro Mirror 3×9 mm – Rectangular Micro Mirror Key Features : High-reflective mirror surface for crystal-clear visibility. Available in four precision mirror sizes. Lightweight ergonomic handle for comfortable control. Durable black anodized aluminum handle. Ideal for indirect vision in confined treatment areas. Autoclavable and reusable. Suitable for endodontic, restorative, and microscopic dental procedures. Applications : Endodontic treatment Root canal procedures Restorative dentistry Intraoral examination Indirect vision under magnification
