



# Prodent Ltd

Quality Dental & Surgical Instruments, Imported from Pakistan, Trusted in Mauritius

---

## Dental Instrument — Dental syringe

Category product list

## Dental syringe

---



### Dental Self-Aspirating Syringe

The Dental Self-Aspirating Syringe is a precision-engineered anesthetic delivery instrument designed for controlled and reliable local anesthesia administration. Manufactured from premium surgical-grade stainless steel with a mirror-polished finish, it offers exceptional durability, corrosion resistance, and comfortable handling for everyday clinical use. Its self-aspirating mechanism enables aspiration without pulling back the thumb ring manually, providing greater control and efficiency during injection procedures.

**Key Features**

- Self-Aspirating Mechanism :** Enables reliable aspiration with minimal hand movement for improved injection control.
- Surgical-Grade Stainless Steel:** Durable, rust-resistant construction for long-lasting performance.
- Ergonomic Finger Rings:** Comfortable grip for enhanced stability and precise handling.
- Mirror-Polished Finish:** Professional appearance with a smooth surface for easy cleaning.
- Reusable & Autoclavable:** Suitable for repeated steam sterilization.
- Universal Cartridge Compatibility** Designed for standard dental anesthetic cartridges and compatible needles.

---



### **Non- Self aspirating syringe**

The Dental Non-Self Aspirating Syringe is a precision-crafted anesthetic delivery instrument designed for the controlled administration of local anesthesia during dental procedures. Manufactured from premium surgical-grade stainless steel with a mirror-polished finish, it delivers exceptional durability, corrosion resistance, and reliable performance for everyday clinical use. Its manual aspiration design provides clinicians with direct control over aspiration through the thumb ring, making it a trusted choice for a wide range of dental treatments.

**Key Features**

- Manual Aspiration Control :** Traditional thumb-ring operation allows precise control during anesthetic administration.
- Surgical-Grade Stainless Steel.** Durable, corrosion-resistant construction for long-term professional use.
- Ergonomic Finger Rings:** Comfortable grip for enhanced stability and accurate handling.
- Mirror-Polished Finish:** Smooth, professional finish for easy cleaning and maintenance.
- Reusable & Autoclavable:** Designed for repeated steam sterilization.
- Standard Compatibility:** Compatible with standard dental anesthetic cartridges and needles.

---