



Prodent Ltd

Quality Dental & Surgical Instruments, Imported from Pakistan, Trusted in Mauritius


Dental Instrument — Elevator

Category product list

Elevator

Curved elevator 7 pcs

Instruments Included : L2S Root Elevator L3S Root Elevator L3C Root Elevator L3IC Root Elevator L3CA Root Elevator L5S Root Elevator L5C Root Elevator Key Benefits : 7-piece professional root elevator set . Multiple working-end patterns for versatile clinical use . Curved designs provide improved access during extraction. Serrated working ends provide enhanced root engagement. Designed for controlled luxation and root elevation. Precision-manufactured working ends. Durable stainless-steel construction . Corrosion-resistant and easy to clean . Suitable for repeated sterilization . Comfortable handling and controlled application . Suitable for professional dental clinics and extraction procedures.



Curved Elevator Kit – 7 Pcs (Micro-Serrated Teething Tips)

The Curved Elevator Kit – 7 Pcs is a premium oral surgery instrument set designed for precise and atraumatic tooth and root extractions. Each elevator features micro-serrated (teething) working tips that provide superior grip on the periodontal ligament and root surface, improving control while reducing slippage during elevation procedures. Manufactured from surgical-grade stainless steel with a mirror-polished finish, this heavy-duty set delivers exceptional durability, corrosion resistance, and long-lasting clinical performance. Supplied in a premium zippered carrying case, the kit keeps all seven elevators organized, protected, and ready for sterilization or transport. Key Features : Complete 7-piece curved elevator set for versatile extraction procedures. Micro-serrated (teething) tips for enhanced grip and precision. Surgical-grade stainless steel construction for superior strength and corrosion resistance. Heavy ergonomic handles for a secure and comfortable grip. Mirror-polished finish for a professional appearance and easy cleaning. Fully autoclavable and reusable. Premium zippered storage case for safe organization and transport. Applications : Atraumatic tooth extractions Root fragment removal Periodontal ligament separation Oral and maxillofacial surgery Routine and surgical extraction procedures




Elevator 13 pieces

13-Piece Root Elevator Kit A complete 13-piece root elevator kit designed for root elevation, loosening, and removal during dental extraction procedures. The kit provides a selection of commonly used root elevators for different clinical requirements. 13 Root Elevators Included : Straight Root Elevator Cryer Root Elevator – Left Cryer Root Elevator – Right Winter Root Elevator – Left Winter Root Elevator – Right Coupland Root Elevator No. 1 Coupland Root Elevator No. 2 Coupland Root Elevator No. 3 Warwick James Root Elevator – Straight Warwick James Root Elevator – Left Warwick James Root Elevator – Right Apical Root Elevator 1 Heidbrink Root Tip Elevator Key product points : 13-piece complete root elevator set. Made from stainless steel Rust-free and corrosion resistant Magnet-free stainless steel Durable construction for repeated clinical use Suitable for dental extraction and root fragment removal Comfortable handling and precise working tips



Luxator Elevator Kit

The 7-Piece Luxator Elevator Kit is designed for controlled periodontal ligament severing, tooth luxation, and atraumatic extraction procedures. The thin, precision-shaped working ends are designed to access the periodontal ligament space and assist in gently loosening teeth before extraction. Instruments Included : Luxator Elevator – Straight Luxator Elevator – Curved Luxator Elevator – Left Luxator Elevator – Right Luxator Elevator – Narrow Luxator Elevator – Medium Luxator Elevator – Wide Key Benefits : 7-piece complete luxator elevator set. Designed for controlled tooth luxation and atraumatic extraction. Thin working ends for improved periodontal ligament access. Different shapes and sizes for versatile clinical applications. Precision-manufactured working tips. Durable stainless-steel construction. Rust-free and corrosion-resistant. Magnet-free stainless steel. Comfortable grip for controlled handling. Suitable for cleaning and repeated sterilization. Ideal for dental clinics and professional extraction procedures



Root elevator 7 pcs

Instruments Included : Straight Root Elevator Cryer Root Elevator – Left Cryer Root Elevator – Right Winter Root Elevator – Left Winter Root Elevator – Right Coupland Root Elevator No. 1 7. Heidbrink Root Tip Elevator Key Benefits : 7-piece complete root elevator set. Designed for root elevation and controlled luxation. Different tip designs for various clinical applications. Durable stainless-steel construction. Corrosion-resistant and easy to clean. Comfortable grip for controlled handling. Suitable for routine dental extraction procedures. Reusable and suitable for sterilization. Professional-quality instruments for dental clinics and practitioners

Straight elevator 5 pcs

Instruments Included : Straight Root Elevator – 2 mm Straight Root Elevator – 3 mm Straight Root Elevator – 4 mm Straight Root Elevator – 5 mm Straight Root Elevator – 6 mm Key Benefits : 5-piece complete straight root elevator set. Different working-end widths for versatile clinical use. Designed for controlled root elevation and luxation. Precision working ends for effective access and handling. Durable stainless-steel construction. Corrosion-resistant and easy to clean. Suitable for repeated sterilization. Comfortable grip for controlled instrument handling



Straight Elevator Kit – 5 Pcs

The Straight Elevator Kit – 5 Pcs is a premium oral surgery instrument set designed for controlled and atraumatic tooth and root extractions. Each elevator features micro-serrated (teething) straight tips that provide superior grip on the periodontal ligament and root surface, enhancing control while minimizing slippage during elevation procedures. Manufactured from surgical-grade stainless steel with a mirror-polished finish, this heavy-duty set delivers exceptional durability, corrosion resistance, and long-lasting clinical performance. Supplied in a premium zippered carrying case, the kit keeps all five elevators securely organized for sterilization, storage, and convenient transport. Key Features : Complete 5-piece straight elevator set for versatile extraction procedures. Micro-serrated (teething) straight tips for enhanced grip and precise control. Surgical-grade stainless steel construction for superior strength and corrosion resistance. Heavy ergonomic handles for a secure and comfortable grip. Mirror-polished finish for easy cleaning and a professional appearance. Fully autoclavable and reusable. Premium zippered storage case for safe organization and transport. Applications : Tooth extractions Root fragment removal Periodontal ligament separation Oral and maxillofacial surgery Routine and surgical extraction procedures
