



Prodent Ltd

Quality Dental & Surgical Instruments, Imported from Pakistan, Trusted in Mauritius

Dental Instrument — Periodontal Instrument

Category product list

Periodontal Instrument

Gracey Curette Set

The 7-Piece Gracey Curette Set is a professional periodontal instrument set designed for scaling, root surface debridement, and periodontal maintenance. Each Gracey curette features a specialized working end designed for specific tooth surfaces and periodontal areas, providing precise and controlled instrumentation. Instruments Included : Gracey Curette 1/2 – Anterior teeth. Gracey Curette 3/4 – Anterior teeth. Gracey Curette 5/6 – Anterior teeth and premolars. Gracey Curette 7/8 – Posterior teeth. Gracey Curette 11/12 – Mesial surfaces of posterior teeth. Gracey Curette 13/14 – Distal surfaces of posterior teeth. Gracey Curette 15/16 – Mesial surfaces of posterior teeth. Key Benefits : 7-piece complete Gracey curette set. Designed for periodontal scaling and root surface debridement. Area-specific working ends for different tooth surfaces. Precision-manufactured for controlled instrumentation. Durable stainless-steel construction. Rust-free and corrosion-resistant. Magnet-free stainless steel. Ergonomic design for comfortable handling. Easy to clean and maintain. Suitable for repeated sterilization. Ideal for dental clinics and periodontal procedures
