



Prodent Ltd

Quality Dental & Surgical Instruments, Imported from Pakistan, Trusted in Mauritius

Dental Rotary & Endodontic Product — Dental Bur

Category product list

Dental Bur



Diamond Burs (Dia Burs)

High-quality diamond burs designed for precise cutting, shaping, preparation, contouring, and finishing during various dental procedures. Available Types of Diamond Burs : Round Bur — For cavity preparation, caries removal, and ccess preparation. Inverted Cone Bur — For cavity preparation, undercuts, and shaping. Straight Fissure Bur — For cutting and shaping tooth struc- ture. Tapered Fissure Bur — For crown preparation, tooth reduction, and shaping. Football / Egg Bur — For occlusal reduction, contouring, and finishing. Flame Bur — For contouring, interproximal finishing, and gingival margin finishing. Needle / Pointed Bur — For detailed finishing, narrow areas, and precise contouring. Cylinder Bur — For flat-surface preparation, reduction, and finishing. Wheel Bur — For occlusal reduction, grooves, and shaping. Chamfer Bur — For crown and bridge preparation and creating chamfer margins. Shoulder Bur — For shoulder-margin crown preparation. End-Cutting Bur — For controlled axial reduction and finishing. Pear-Shaped Bur — For cavity preparation and restorative procedures. Lentulo / Special-Shaped Burs — Available for specific dental laboratory and clinical applications depending on the design. Finishing & Polishing Burs — For smooth finishing and refinement of restorations. Key Benefits : Precision cutting and preparation. Smooth and controlled performance. Suitable for restorative and prosthetic procedures. Multiple shapes and sizes available. Diamond-coated working surface for efficient cutting. Compatible with standard high-speed dental handpieces. Durable construction for reliable clinical use