



Prodent Ltd

Quality Dental & Surgical Instruments, Imported from Pakistan, Trusted in Mauritius

Orthodontic

Full catalog

Orthodontic instrument



Band Pusher / Band Seater

Designed for seating and adapting orthodontic bands during placement.



Bracket Holder

Used for securely holding and positioning brackets during placement. Stainless steel. Rust-free. Magnet-free.



Bracket Positioning Instrument

Designed to assist with accurate positioning and placement of orthodontic brackets. Stainless steel. Rust-free. Magnet-free.



Ligature Director

Designed for guiding and positioning ligatures around orthodontic brackets. Stainless steel. Rust-free. Magnet-free.



Ligature Tucking Instrument

Used for tucking and securing ligature wires around orthodontic brackets. Stainless steel. Rust-free. Magnet-free.



Orthodontic Explorer / Probe

Used for examination and access during orthodontic procedures. Stainless steel. Rust-free. Magnet-free.



Orthodontic Scaler

Designed for cleaning and removing deposits around orthodontic appliances. Stainless steel Rust-free Magnet-free.

Orthodontic kit



21- piece ortho kit

21-Piece Complete Orthodontic Instrument Kit A complete 21-piece orthodontic instrument kit designed to provide dental professionals with a convenient selection of essential instruments for orthodontic procedures. The instruments are made from stainless steel, rust-free, and magnet-free ,making this kit a practical choice for professional dental and orthodontic use. Key Features : Complete 21-piece orthodontic kit. Stainless steel. Rust-free. Magnet-free. Essential instruments for orthodontic procedures. Suitable for professional dental and orthodontic clinics.

Pliers



Adams Plier

Designed for bending and forming orthodontic wires, particularly for removable appliances. Stainless steel. Rust-free. Magnet-free.



Archwire Bending Plier

Designed for making precise bends and adjustments to orthodontic archwires. Stainless steel. Rust-free. Magnet-free.



Band Removing Plier

Used for efficient removal of orthodontic bands. Stainless steel. Rust-free. Magnet-free.



Bird Beak Plier

Used for bending and forming orthodontic wires into precise shapes and loops. Stainless steel. Rust-free. Magnet-free.



Bracket Removing Plier

Designed for controlled removal of orthodontic brackets. Stainless steel, rust-free, and magnet-free

Distal End Cutter Plier

Used for cutting the distal end of orthodontic arch wires after placement and adjustment. Stainless steel Rust-free Magnet-free



How Plier

Used for gripping, placing, and manipulating orthodontic wires and components. Stainless steel. Rust-free. Magnet-free.



Ligature Cutter Plier

Designed for cutting ligature wires efficiently during orthodontic procedures. Stainless steel. Rust-free. Magnet-free.



Ligature Director / Cutter

Used for guiding, adjusting, and managing ligature wires during orthodontic treatment. Stainless steel Rust-free Magnet-free.



Mathieu Orthodontic Plier

Used for handling and placing ligatures, elastics, and other small orthodontic components. Stainless steel. Rust-free. Magnet-free.

Three-Prong Plier

Designed for controlled bending and adjustment of orthodontic wires. Stainless steel. Rust-free. Magnet-free.



Tweed Arch Forming Plier

Used for shaping and forming orthodontic archwires to the required arch form. Stainless steel. Rust-free. Magnet-free.



Weingart Utility Plier

Designed for gripping, guiding, placing, and removing arch wires during orthodontic procedures. Stainless steel Rust-free Magnet-free
