



Prodent Ltd

Quality Dental & Surgical Instruments, Imported from Pakistan, Trusted in Mauritius

Orthodontic — Orhtodontic instrument

Category product list

Orhtodontic instrument



Band Pusher / Band Seater

Designed for seating and adapting orthodontic bands during placement.



Bracket Holder

Used for securely holding and positioning brackets during placement. Stainless steel. Rust-free. Magnet-free.



Bracket Positioning Instrument

Designed to assist with accurate positioning and placement of orthodontic brackets. Stainless steel. Rust-free. Magnet-free.



Ligature Director

Designed for guiding and positioning ligatures around orthodontic brackets. Stainless steel. Rust-free. Magnet-free.



Ligature Tucking Instrument

Used for tucking and securing ligature wires around orthodontic brackets. Stainless steel. Rust-free. Magnet-free.



Orthodontic Explorer / Probe

Used for examination and access during orthodontic procedures. Stainless steel. Rust-free. Magnet-free.



Orthodontic Scaler

Designed for cleaning and removing deposits around orthodontic appliances. Stainless steel Rust-free Magnet-free.
