



Prodent Ltd

Quality Dental & Surgical Instruments, Imported from Pakistan, Trusted in Mauritius

Orthodontic — Pliers

Category product list

Pliers



Adams Plier

Designed for bending and forming orthodontic wires, particularly for removable appliances. Stainless steel. Rust-free. Magnet-free.



Archwire Bending Plier

Designed for making precise bends and adjustments to orthodontic archwires. Stainless steel. Rust-free. Magnet-free.



Band Removing Plier

Used for efficient removal of orthodontic bands. Stainless steel. Rust-free. Magnet-free.



Bird Beak Plier

Used for bending and forming orthodontic wires into precise shapes and loops. Stainless steel. Rust-free. Magnet-free.



Bracket Removing Plier

Designed for controlled removal of orthodontic brackets. Stainless steel, rust-free, and magnet-free

Distal End Cutter Plier

Used for cutting the distal end of orthodontic arch wires after placement and adjustment. Stainless steel Rust-free Magnet-free



How Plier

Used for gripping, placing, and manipulating orthodontic wires and components. Stainless steel. Rust-free. Magnet-free.



Ligature Cutter Plier

Designed for cutting ligature wires efficiently during orthodontic procedures. Stainless steel. Rust-free. Magnet-free.



Ligature Director / Cutter

Used for guiding, adjusting, and managing ligature wires during orthodontic treatment. Stainless steel Rust-free Magnet-free.



Mathieu Orthodontic Plier

Used for handling and placing ligatures, elastics, and other small orthodontic components. Stainless steel. Rust-free. Magnet-free.

Three-Prong Plier

Designed for controlled bending and adjustment of orthodontic wires. Stainless steel. Rust-free. Magnet-free.



Tweed Arch Forming Plier

Used for shaping and forming orthodontic archwires to the required arch form. Stainless steel. Rust-free. Magnet-free.



Weingart Utility Plier

Designed for gripping, guiding, placing, and removing arch wires during orthodontic procedures. Stainless steel Rust-free Magnet-free
